

# Certificate of Analysis for HM-277D

## Genomic DNA from Microbial Mock Community B (Staggered, High Concentration), v5.2H, for Whole Genome Shotgun Sequencing

Catalog No. HM-277D

**Product Description:** A mixture of genomic DNA from 20 bacterial strains containing staggered ribosomal RNA operon counts (10,000 – 10,000,000 operons per organism per  $\mu$ L). **Note: The label for HM-277D is incorrect. HM-277D contains genomic DNA from microbial mock community B and not microbial mock community A.**

Lot<sup>1,2</sup>: 60257284

Manufacturing Date: 31AUG2011

TEST	SPECIFICATIONS	RESULTS
DNA Sequencing of Individual 16S Ribosomal RNA Genes from Mock Community B (~ 1500 base pairs)	<p>Consistent with <i>Acinetobacter baumannii</i></p> <p>Consistent with <i>Actinomyces odontolyticus</i></p> <p>Consistent with <i>Bacillus cereus</i></p> <p>Consistent with <i>Bacteroides vulgatus</i></p> <p>Consistent with <i>Clostridium beijerinckii</i></p> <p>Consistent with <i>Deinococcus radiodurans</i></p> <p>Consistent with <i>Enterococcus faecalis</i></p> <p>Consistent with <i>Escherichia coli</i></p> <p>Consistent with <i>Helicobacter pylori</i></p> <p>Consistent with <i>Lactobacillus gasseri</i></p> <p>Consistent with <i>Listeria monocytogenes</i></p> <p>Consistent with <i>Neisseria meningitidis</i></p> <p>Consistent with <i>Propionibacterium acnes</i></p> <p>Consistent with <i>Pseudomonas aeruginosa</i></p> <p>Consistent with <i>Rhodobacter sphaeroides</i></p> <p>Consistent with <i>Staphylococcus aureus</i></p> <p>Consistent with <i>Staphylococcus epidermidis</i></p> <p>Consistent with <i>Streptococcus agalactiae</i></p> <p>Consistent with <i>Streptococcus mutans</i></p> <p>Consistent with <i>Streptococcus pneumoniae</i></p>	<p>Consistent with <i>Acinetobacter baumannii</i><sup>§</sup></p> <p>Consistent with <i>Actinomyces odontolyticus</i><sup>§</sup></p> <p>Consistent with <i>Bacillus cereus</i><sup>§</sup></p> <p>Consistent with <i>Bacteroides vulgatus</i><sup>§</sup></p> <p>Consistent with <i>Clostridium beijerinckii</i><sup>†</sup></p> <p>Consistent with <i>Deinococcus radiodurans</i><sup>§</sup></p> <p>Consistent with <i>Enterococcus faecalis</i><sup>§</sup></p> <p>Consistent with <i>Escherichia coli</i><sup>†</sup></p> <p>Consistent with <i>Helicobacter pylori</i><sup>†</sup></p> <p>Consistent with <i>Lactobacillus gasseri</i><sup>†</sup></p> <p>Consistent with <i>Listeria monocytogenes</i><sup>§</sup></p> <p>Consistent with <i>Neisseria meningitidis</i><sup>†</sup></p> <p>Consistent with <i>Propionibacterium acnes</i><sup>§,3</sup></p> <p>Consistent with <i>Pseudomonas aeruginosa</i><sup>£</sup></p> <p>Consistent with <i>Rhodobacter sphaeroides</i><sup>£</sup></p> <p>Consistent with <i>Staphylococcus aureus</i><sup>§,4</sup></p> <p>Consistent with <i>Staphylococcus epidermidis</i><sup>§,4</sup></p> <p>Consistent with <i>Streptococcus agalactiae</i><sup>§</sup></p> <p>Consistent with <i>Streptococcus mutans</i><sup>§</sup></p> <p>Consistent with <i>Streptococcus pneumoniae</i><sup>§</sup></p>
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)

# Certificate of Analysis for HM-277D

TEST	SPECIFICATIONS	RESULTS
<b>Theoretical DNA Concentration for Individual Members of Mock Community B [based on number of ribosomal RNA (rRNA) operons input DNA]</b> <i>Acinetobacter baumannii</i> - 100,000 operons <i>Actinomyces odontolyticus</i> - 10,000 operons <i>Bacillus cereus</i> - 1,000,000 operons <i>Bacteroides vulgatus</i> - 10,000 operons <i>Clostridium beijerinckii</i> - 1,000,000 operons <i>Deinococcus radiodurans</i> - 10,000 operons <i>Enterococcus faecalis</i> - 10,000 operons <i>Escherichia coli</i> - 10,000,000 operons <i>Helicobacter pylori</i> - 100,000 operons <i>Lactobacillus gasseri</i> - 100,000 operons <i>Listeria monocytogenes</i> - 100,000 operons <i>Neisseria meningitidis</i> - 100,000 operons <i>Propionibacterium acnes</i> - 100,000 operons <i>Pseudomonas aeruginosa</i> - 1,000,000 operons <i>Rhodobacter sphaeroides</i> - 10,000,000 operons <i>Staphylococcus aureus</i> - 1,000,000 operons <i>Staphylococcus epidermidis</i> - 10,000,000 operons <i>Streptococcus agalactiae</i> - 1,000,000 operons <i>Streptococcus mutans</i> - 10,000,000 operons <i>Streptococcus pneumoniae</i> - 10,000 operons	Report results	82 pg/μL <i>Acinetobacter baumannii</i> <sup>§</sup> 10 pg/μL <i>Actinomyces odontolyticus</i> <sup>§</sup> 450 pg/μL <i>Bacillus cereus</i> <sup>§</sup> 7.6 pg/μL <i>Bacteroides vulgatus</i> <sup>§</sup> 440 pg/μL <i>Clostridium beijerinckii</i> <sup>‡</sup> 10 pg/μL <i>Deinococcus radiodurans</i> <sup>§</sup> 7.0 pg/μL <i>Enterococcus faecalis</i> <sup>§</sup> 6.8 ng/μL <i>Escherichia coli</i> <sup>£</sup> 86 pg/μL <i>Helicobacter pylori</i> <sup>†</sup> 32 pg/μL <i>Lactobacillus gasseri</i> <sup>†</sup> 50 pg/μL <i>Listeria monocytogenes</i> <sup>§</sup> 58 pg/μL <i>Neisseria meningitidis</i> <sup>†</sup> 88 pg/μL <i>Propionibacterium acnes</i> <sup>§</sup> 1.6 ng/μL <i>Pseudomonas aeruginosa</i> <sup>£</sup>  14 ng/μL <i>Rhodobacter sphaeroides</i> <sup>£</sup>  590 pg/μL <i>Staphylococcus aureus</i> <sup>§</sup> 5.1 ng/μL <i>Staphylococcus epidermidis</i> <sup>§</sup>  32 pg/μL <i>Streptococcus agalactiae</i> <sup>§</sup>  4.1 ng/μL <i>Streptococcus mutans</i> <sup>§</sup> 5.5 pg/μL <i>Streptococcus pneumoniae</i> <sup>§</sup>
<b>Total Amount of DNA per vial</b>	≥ 50 ng per μL	56 ng per μL
<b>Functional Activity by PCR Amplification</b> 16S ribosomal RNA gene	~ 1500 base pair amplicon	~ 1500 base pair amplicon (Figure 1)
<b>OD<sub>260</sub>/OD<sub>280</sub> Ratio</b>	Report results	1.9
<b>Bacterial Inactivation</b> 10% of total yield plated on Tryptic Soy agar with 5% sheep blood <sup>5</sup>	No viable bacteria detected	No viable bacteria detected

<sup>1</sup>Preparation and QC testing (with the exception of Bacterial Inactivation) were performed by Baylor College of Medicine in Houston, Texas.


<sup>2</sup>Genomic DNA was extracted using the following methods: <sup>§</sup>SDS Lysis, CsCl, <sup>‡</sup>Modified SDS Lysis, CsCl, <sup>£</sup>Triton Lysis, CsCl and <sup>†</sup>Omega E.Z.N.A.<sup>®</sup> Bacterial DNA Kit.

<sup>3</sup>Also consistent with other *Propionibacterium* species

<sup>4</sup>Also consistent with other *Staphylococcus* species

<sup>5</sup>7 days at 37°C under both anaerobic atmosphere (80% N<sub>2</sub>:10% CO<sub>2</sub>:10% H<sub>2</sub>) and aerobic atmospheric conditions

# Certificate of Analysis for HM-277D



Lane 1: 1 Kb DNA Ladder (Fermentas, Cat. No. SM0311)

Lane 2: 200 ng of gDNA HM-277D

Lane 3: PCR of 16S ribosomal RNA gene from HM-277D

**Date:** 07 OCT 2014

**Signature:**

**Title:**

Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by ATCC® and the contractor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

